

Intermittent Catheters and Trays

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician

INDICATIONS FOR USE: For urological use only.

WARNINGS

- On latex and red rubber catheter, do not use ointments or lubricants having a petrolatum base. They will damage latex.
- Visually inspect the product for any imperfections or surface deterioration prior to use.
- After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.
- Some Intermittent Catheters may contain natural rubber latex which may cause allergic reaction. Please check package.
- Contains or presence of Phthalates: Di(2-ethylhexyl)phthalate (DEHP) is a plasticizer used in some polyvinyl chloride medical devices. DEHP has been shown to produce a range of adverse effects in experimental animals, notably liver toxicity and testicular atrophy. Although the toxic and carcinogenic effects of DEHP have been well established in experimental animals, the ability of this compound to produce adverse effects in humans is controversial. There is no evidence that neonates, infants, pregnant and breast feeding women exposed to DEHP experience any related adverse effects. However, a lack of evidence of causation between DEHP-PVC and any disease or adverse effect does not mean that there are no risks.

PRECAUTIONS(BARD® TOUCHLESS® Intermittent Catheter Only)

- It is recommended that the collection bag be held. The catheter could possibly separate from the collection bag when urine increases weight of the bag.

NOTES

- This is a single use device.
- Do not re-sterilize any portion of this device.
- Reuse and/or repackaging may create a risk of patient or user infection, compromise the structural integrity and/or essential material and design characteristics of device, which may lead to device failure, and/or lead to injury, illness or death of the patient.
- Single patient use, do not reuse.
- Sterile unless package is opened or damaged.
- Storage: Store catheters at room temperature away from direct exposure to light, preferably in the original box.

For the latest information, always check the "Instructions for Use" that comes packaged with the product. If you are a consumer seeking more information, please consult your healthcare provider.