

COMFASURE® Catheter Retainer Strap

BARD® Comprehensive Care

Supporting your Collection System

Extra security for total confidence.

COMFASURE® fits comfortably to the leg or abdomen when wearing a catheter, to give you the security, confidence and reassurance you need to enjoy the fullest possible life.


Easy to Use

Softly elasticated COMFASURE® has an anti-slip backing. It's unique, easy-to-operate retainer clip secures with a 'click' - supporting your Foley catheter or drainage bag tubing neatly without crushing. You have the assurance to move freely, knowing your catheter is safely held in position.

Easy to Choose

COMFASURE® is available in three sizes. To determine the correct size needed, measure the circumference of the widest point of the thigh (A) or abdomen.

COMFASURE® is available on prescription in packs of five and is available on prescription from your Nurse or Doctor.



Fitting Guide

COMFASURE® Size	Circumference	Product Code
Adult	Up to 65cm	AD3403
Small	Up to 35cm	SM3404
Abdominal#		AB3405



FREEPHONE 0800 0121 699

Fax: 01903 875 085

Email: ukcustomer care@crbard.com

Web: www.scripteasy.co.uk

Brightwake Limited, Lowmoor Business Park, Kirkby-in-Ashfield, Nottinghamshire, NG17 7JZ

Distributed by Bard Limited, Forest House, Tilgate Forest Business Park, Brighton Road, Crawley, West Sussex, RH11 9BP.
Tel: +44 (0) 1293 527 888 Fax: +44 (0) 1293 552 428

Bard and Comfasure are registered trademarks of C. R. Bard, Inc., or an affiliate. © 2012, C. R. Bard, Inc. All rights reserved. Consult product labels and inserts for any indications, contraindications, hazards, warnings cautions and instructions for use. #Abdominal COMFASURE® can be cut and adjusted to fit a larger leg. 0213/4082

NURSE ADVICE LINE FREEPHONE 0808 168 4048

www.bardmedical.co.uk

BARD | MEDICAL