

Magic³ Touchless™ New Features

HYDROPHILIC INTERMITTENT CATHETER SYSTEM

	MAGIC ³ Closed System	MAGIC ³ Touchless™ Closed System
The Catheter		
MAGIC ³ all-silicone catheter	•	•
Unique three-layer design	•	•
Water-activated hydrophilic coating	•	•
Does not contain latex, DEHP, or phthalates	•	•
Tapered tip	•	•
Four comfort-sized drainage eyes	•	•
No-touch insertion method	•	•
French sizes	8 - 18	8 - 16
The Bag		
1100cc collection bag	•	•
Contains sterile water packet	•	•
Built-in handle	•	•
Instruction illustrations directly on bag	•	•
Tear feature for drainage	•	•
Tapered neck bag design	•	•
Additional Features		
Designed to help reduce the risk of UTIs	•	•
PVI or BZK accessory kit options	•	•
Pre-lubricated introducer tip	•	•
SURE-GRIP™ sleeve	•	•
Catheter cap ring	Clear	Opaque (White)
Quantity	50/case	30/case

The PERSONAL CATHETER® Closed System includes all of the features of the MAGIC³ Closed System, excluding the three-layer silicone hydrophilic catheter. The PERSONAL CATHETER® contains a single-layer all-silicone hydrophilic catheter.



Care. Comfort. Choice.™

Bard Medical Division
Covington, GA
1.800.526.4455
www.bardcare.com

Summary of Indications & Warnings:

Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.

Bard, Magic³, Magic³ Touchless, Personal Catheter and Sure-Grip are trademarks and/or registered trademarks of C.R. Bard Inc. ©2016 C.R. Bard Inc. All Rights Reserved. Printed in U.S.A.
1603-31 R04/16 BMD/TCHL/0416/0005