

Everything is  
taken care of.

**Spirit**<sup>®</sup>

LEG BAGS & BED BAGS

**BARD** | **CARE**

Care. Comfort. Choice.

**NEW**

# Spirit®

LEG BAGS & BED BAGS

MORE THAN  
**20%**  
LESS EXPENSIVE  
THAN UK MARKET  
LEADERS\*

## Available on Drug Tariff

- 1 Ridge Connector**  
Providing a secure connection
- 2 Latex Free Sterile Gloves**  
Protecting against cross-infection
- 3 Needle-free Sample Port**  
Facilitating safe sampling
- 4 180° Lever Tap**  
Easy to operate security
- 5 Leg Straps**  
Offering added security
- 6 Docking Station**  
Minimising infection risks via Bed Bag outlet tube



### SPIRIT® Leg Bags

### Ordering Information

CODE	VOLUME	TUBE TYPE	TUBE LENGTH	QTY	PRICE
DSL350S	350mL	Short Tube	10cm	10	£24.00
DSL500S	500mL	Short Tube	10cm	10	£24.00
DSL750S	750mL	Short Tube	10cm	10	£24.00
DSL350L	350mL	Long Tube	30cm	10	£24.00
DSL500L	500mL	Long Tube	30cm	10	£24.00
DSL750L	750mL	Long Tube	30cm	10	£24.00
DSL750V	750mL	Variable Tube	60cm	10	£24.00

### SPIRIT® Bed Bags

CODE	VOLUME	TUBE TYPE	TUBE LENGTH	QTY	PRICE
DSBB2000	2000mL	Long Tube	120cm	10	£11.00



\*Compared to pricing in January 2017 National Health Service England and Wales Drug Tariff ISBN: 978-0-11-784296-0



**FREEPHONE 0800 470 1930**  
HEALTHCARE PROFESSIONALS

**NURSE ADVICE LINE**  
**FREEPHONE**

Fax  
**01903 875 085**

Email  
**ukcustomercare@crbard.com**

Web  
**www.scripteasy.co.uk**

**0808 168 4048**

**Bard Limited**, Forest House, Tilgate Forest Business Park, Brighton Road, Crawley, West Sussex, RH11 9BP, UK  
Telephone: **01293 527 888** Fax: **01293 552 428** Bard Customer Care: **01293 529 555** Web: **www.bardcare.uk**



Please consult product labels and insert for any indications, contraindications, hazards, warnings, cautions and directions for use. Bard, Script-easy and Spirit are trademarks and/or registered trademarks of C. R. Bard, Inc. © 2016 C. R. Bard, Inc. All Rights Reserved. Script-easy® is the prescription dispensing service operated by Bard Limited. All Rights Reserved. 1116/5006 17-017